

COVID-19 Rapid Antigen Testing



Overview

- Abbott BINAXNOW
- Anterior nares swab
- Detects viral antigen
- Results available within 15 minutes
- 600,000 kits may be allocated to Fire/EMS by Testing Task Force




VIDEO: <https://www.globalpointofcare.abbott/en/product-details/navica-binaxnow-covid-19-us.html#>

Overview

- Covered under MIEMSS CLIA waiver
- No physician orders required
- Requires reporting of all results via CRISP



Antigen Testing vs. PCR





MJ Mina et al. N Engl J Med 2020. DOI: 10.1056/NEJMp2025631



The NEW ENGLAND
JOURNAL of MEDICINE

Test Results

- Symptomatic patients within 7 days of onset
 - Sensitivity 97%
 - Specificity 98.5%
- Invalid result – repeat antigen test



Use Cases

- Evaluate a symptomatic person**
- Contact tracing / investigation / surveillance
- Screening (with frequency)

Follow-up PCR Testing

	Symptomatic person	Asymptomatic person
Antigen test POS	Positive	*Obtain PCR test
Antigen test NEG	*Serial Testing	Negative

Truth

	<i>Disease (number)</i>	<i>Non Disease (number)</i>	<i>Total (number)</i>
<i>Positive (number)</i>	A <i>(True Positive)</i>	B <i>(False Positive)</i>	$T_{\text{Test Positive}}$
<i>Negative (number)</i>	C <i>(False Negative)</i>	D <i>(True Negative)</i>	$T_{\text{Test Negative}}$
	T_{Disease}	$T_{\text{Non Disease}}$	Total

Specificity = True Negative / (True Negative + False Positive)

Positive Predictive Value = True Positive / (True Positive + False Positive)

Positive Predictive Value** if disease prevalence among the tested population is **10%**

	Has COVID-19	Doesn't Have COVID
Test Result Positive	9.7	1.8
Test Result Negative	0.3	88.2
Total (if 100 tested)	10	90

** Assuming test sensitivity = 97%, and specificity = 98%

Positive Predictive Value = 84%

(84% of people with a positive result will actually have COVID-19)

Positive Predictive Value** if disease prevalence among the tested population is 5%

	Has COVID-19	Doesn't Have COVID
Test Result Positive	4.85	1.9
Test Result Negative	.15	93.1
Total (if 100 tested)	5	95

** Assuming test sensitivity = 97%, and specificity = 98%

Positive Predictive Value = 72%

(72% of people with a positive result will actually COVID-19)

Currently, across the state, PCR testing of EMS/Fire personnel is Positive 8-9%. Presumably, these people are being tested because of suspicion (i.e., higher pre-test probability) or concern as a contact. Lower prevalence, and decreased positive predictive value, should be anticipated if testing was to be used for screening. That's not bad, but something to plan for.

Process

- EMSOP submits a written plan and is approved by MIEMSS
 - Process for manually reporting all COVID-19 antigen tests (pos and neg) into CRISP
 - Process for ensuring clinicians who require a PCR test receive one within 24-48 hours of positive antigen result
 - Process to ensure weekly COVID-19 antigen inventory is submitted through SmartSheet form

Logistics

- Test kits distributed in sets of 40
- Same process as currently used for PCR test kits
- Six month expiration (rolling) from date of arrival at MDH

Training

- Abbott offers virtual super-user training
- Training resources

<https://www.globalpointofcare.abbott/en/support/product-installation-training/navica-brand/navica-binaxnow-ag-training.html>